

Peptidyl-prolyl cis-trans isomerase B

Catalog No: #AP86950



Package Size: #AP86950-1 50ug #AP86950-2 100ug #AP86950-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peptidyl-prolyl cis-trans isomerase B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PPlase B,CYP-S1,Cyclophilin B,Rotamase B,S-cyclophilin,SCYLP
Accession No.	Uniprot ID:P24369
Uniprot	P24369
GeneID	19035;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ppib

Note: This product is for in vitro research use only