

# Serine/threonine-protein phosphatase PP1-alpha catalytic subunit

Catalog No: #AP86956

Package Size: #AP86956-1 50ug #AP86956-2 100ug #AP86956-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine/threonine-protein phosphatase PP1-alpha catalytic subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PP-1A
Accession No.	Uniprot ID:Q3T0E7
Uniprot	Q3T0E7
GeneID	516175;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PPP1CA

Note: This product is for in vitro research use only