

Serum paraoxonase/arylesterase 2

Catalog No: #AP86958



Package Size: #AP86958-1 50ug #AP86958-2 100ug #AP86958-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serum paraoxonase/arylesterase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	PON 2,Aromatic esterase 2,A-esterase 2,Serum aryldialkylphosphatase 2
Accession No.	Uniprot ID:Q90952
Uniprot	Q90952
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PON2

Note: This product is for in vitro research use only