

Peroxisome proliferator-activated receptor delta

Catalog No: #AP86960

Package Size: #AP86960-1 50ug #AP86960-2 100ug #AP86960-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peroxisome proliferator-activated receptor delta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	PPAR-delta,Nuclear receptor subfamily 1 group C member 2,Peroxisome proliferator-activated receptor beta,PPAR-beta,NR1C2,PPARB
Accession No.	Uniprot ID:Q0ZQA8
Uniprot	Q0ZQA8
GeneID	481756;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PPARD

Note: This product is for in vitro research use only