

Alkaline phosphatase, tissue-nonspecific isozyme

Catalog No: #AP86961



Package Size: #AP86961-1 50ug #AP86961-2 100ug #AP86961-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alkaline phosphatase, tissue-nonspecific isozyme
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AP-TNAP,Alkaline phosphatase liver/bone/kidney isozyme
Accession No.	Uniprot ID:P08289
Uniprot	P08289
GeneID	25586;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Alpl

Note: This product is for in vitro research use only