

Prolactin-inducible protein homolog

Catalog No: #AP86963

Package Size: #AP86963-1 50ug #AP86963-2 100ug #AP86963-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prolactin-inducible protein homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Prolactin-induced protein
Accession No.	Uniprot ID:P60986
Uniprot	P60986
GeneID	613853;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PIP

Note: This product is for in vitro research use only