

Phospholipid scramblase 1

Catalog No: #AP86964



Package Size: #AP86964-1 50ug #AP86964-2 100ug #AP86964-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospholipid scramblase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PL scramblase 1,Ca(2+)-dependent phospholipid scramblase 1,Transplantability-associated protein 1,NOR1,Tra1b,Tras1
Accession No.	Uniprot ID:Q9JJ00
Uniprot	Q9JJ00
GeneID	22038;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Plscr1

Note: This product is for in vitro research use only