

Polycystin-1

Catalog No: #AP86965



Package Size: #AP86965-1 50ug #AP86965-2 100ug #AP86965-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Polycystin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PC1,Autosomal dominant polycystic kidney disease 1 protein
Accession No.	Uniprot ID:P98161
Uniprot	P98161
GeneID	5310;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PKD1

Note: This product is for in vitro research use only