

# cAMP and cAMP-inhibited cGMP 3',5'-cyclic phosphodiesterase 10A



Catalog No: #AP86966

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP86966-1 50ug #AP86966-2 100ug #AP86966-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	cAMP and cAMP-inhibited cGMP 3',5'-cyclic phosphodiesterase 10A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:Q9Y233
Uniprot	Q9Y233
GeneID	10846;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: PDE10A

Note: This product is for in vitro research use only