

Hemojuvelin

Catalog No: #AP86968



Package Size: #AP86968-1 50ug #AP86968-2 100ug #AP86968-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hemojuvelin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Hemochromatosis type 2 protein homolog,RGM domain family member C,Rgmc
Accession No.	Uniprot ID:Q8N7M5
Uniprot	Q8N7M5
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hfe2

Note: This product is for in vitro research use only