

Putative testis-specific prion protein

Catalog No: #AP86969



Package Size: #AP86969-1 50ug #AP86969-2 100ug #AP86969-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Putative testis-specific prion protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein M8
Accession No.	Uniprot ID:Q86SH4
Uniprot	Q86SH4
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PRNT

Note: This product is for in vitro research use only