

Proteasome subunit alpha type-1

Catalog No: #AP86970



Package Size: #AP86970-1 50ug #AP86970-2 100ug #AP86970-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proteasome subunit alpha type-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	30 kDa prosomal protein,PROS-30,Macropain subunit C2,Multicatalytic endopeptidase complex subunit C2,Proteasome component C2,Proteasome nu chain,HC2,NU,PROS30,PSC2
Accession No.	Uniprot ID:P25786
Uniprot	P25786
GeneID	5682;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PSMA1

Note: This product is for in vitro research use only