

Proteasome subunit alpha type-1

Catalog No: #AP86973



Package Size: #AP86973-1 50ug #AP86973-2 100ug #AP86973-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proteasome subunit alpha type-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Macropain subunit C2, Multicatalytic endopeptidase complex subunit C2, Proteasome component C2, Proteasome nu chain
Accession No.	Uniprot ID:Q9R1P4
Uniprot	Q9R1P4
GeneID	26440;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PsmA1

Note: This product is for in vitro research use only