

Advanced glycosylation end product-specific receptor

Catalog No: #AP86975



Package Size: #AP86975-1 50ug #AP86975-2 100ug #AP86975-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Advanced glycosylation end product-specific receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Receptor for advanced glycosylation end products,RAGE
Accession No.	Uniprot ID:Q28173
Uniprot	Q28173
GeneID	280986;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AGER

Note: This product is for in vitro research use only