

Pituitary adenylate cyclase-activating polypeptide

Catalog No: #AP86976



Package Size: #AP86976-1 50ug #AP86976-2 100ug #AP86976-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pituitary adenylate cyclase-activating polypeptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PACAP,PACAP
Accession No.	Uniprot ID:Q29W19
Uniprot	Q29W19
GeneID	615187;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADCYAP1

Note: This product is for in vitro research use only