

Renalase

Catalog No: #AP86981



Package Size: #AP86981-1 50ug #AP86981-2 100ug #AP86981-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Renalase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Monoamine oxidase-C, MAO-C, alpha-NAD(P)H oxidase/anomerase
Accession No.	Uniprot ID:Q5U2W9
Uniprot	Q5U2W9
GeneID	361751;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Rnls

Note: This product is for in vitro research use only