

Polycystin-2

Catalog No: #AP86982



Package Size: #AP86982-1 50ug #AP86982-2 100ug #AP86982-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Polycystin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Polycystic kidney disease 2 protein homolog, Transient receptor potential cation channel subfamily P member 2, TRPP2
Accession No.	Uniprot ID:O35245
Uniprot	O35245
GeneID	18764;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pkd2

Note: This product is for in vitro research use only