

## Peroxisome proliferator-activated receptor alpha

Catalog No: #AP86983



Package Size: #AP86983-1 50ug #AP86983-2 100ug #AP86983-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Peroxisome proliferator-activated receptor alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PPAR-alpha,Nuclear receptor subfamily 1 group C member 1,Nr1c1,Ppar
Accession No.	Uniprot ID:P37230
Uniprot	P37230
GeneID	25747;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ppara

Note: This product is for in vitro research use only