

# Retinoic acid receptor gamma

Catalog No: #AP86987



Package Size: #AP86987-1 50ug #AP86987-2 100ug #AP86987-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Retinoic acid receptor gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	RAR-gamma,Nuclear receptor subfamily 1 group B member 3,Nr1b3
Accession No.	Uniprot ID:P18911
Uniprot	P18911
GeneID	19411;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Rarg

Note: This product is for in vitro research use only