

Thioredoxin-dependent peroxide reductase, mitochondrial

Catalog No: #AP86988

Package Size: #AP86988-1 50ug #AP86988-2 100ug #AP86988-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thioredoxin-dependent peroxide reductase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PRx III, Peroxiredoxin-3, PRX-3
Accession No.	Uniprot ID:Q9Z0V6
Uniprot	Q9Z0V6
GeneID	64371;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Prdx3

Note: This product is for in vitro research use only