

## 26S protease regulatory subunit 6A

Catalog No: #AP86989



Package Size: #AP86989-1 50ug #AP86989-2 100ug #AP86989-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	26S protease regulatory subunit 6A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	26S proteasome AAA-ATPase subunit RPT5, Proteasome 26S subunit ATPase 3, Tat-binding protein 1, TBP-1, Tbp1
Accession No.	Uniprot ID:O88685
Uniprot	O88685
GeneID	19182;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Psmc3

Note: This product is for in vitro research use only