

Oxysterols receptor LXR-beta

Catalog No: #AP86990



Package Size: #AP86990-1 50ug #AP86990-2 100ug #AP86990-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oxysterols receptor LXR-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Liver X receptor beta,Nuclear receptor subfamily 1 group H member 2,Retinoid X receptor-interacting protein No.15,Ubiquitously-expressed nuclear receptor,Lxrb,Rip15,Unr,Unr2
Accession No.	Uniprot ID:Q60644
Uniprot	Q60644
GeneID	22260;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nr1h2

Note: This product is for in vitro research use only