

Neurexin-1

Catalog No: #AP86992



Package Size: #AP86992-1 50ug #AP86992-2 100ug #AP86992-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neurexin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Neurexin I-alpha,Neurexin-1-alpha,Kiaa0578
Accession No.	Uniprot ID:Q9CS84
Uniprot	P0DI97
GeneID	18189;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nrxn1

Note: This product is for in vitro research use only