

Phospholipase A2

Catalog No: #AP86995



Package Size: #AP86995-1 50ug #AP86995-2 100ug #AP86995-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospholipase A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Group IB phospholipase A2, Phosphatidylcholine 2-acylhydrolase 1B
Accession No.	Uniprot ID:P00593
Uniprot	P00593
GeneID	282457;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PLA2G1B

Note: This product is for in vitro research use only