

Protein phosphatase methylesterase 1

Catalog No: #AP86998



Package Size: #AP86998-1 50ug #AP86998-2 100ug #AP86998-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein phosphatase methylesterase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PME-1,PME1,PP2593,PRO0750
Accession No.	Uniprot ID:Q9Y570
Uniprot	Q9Y570
GeneID	51400;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PPME1

Note: This product is for in vitro research use only