

Peroxiredoxin-5, mitochondrial

Catalog No: #AP86999



Package Size: #AP86999-1 50ug #AP86999-2 100ug #AP86999-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peroxiredoxin-5, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Antioxidant enzyme B166,AOEB166,Liver tissue 2D-page spot 2D-0014IV,PLP,Peroxiredoxin V,Prx-V,Peroxisomal antioxidant enzyme,Thioredoxin peroxidase PMP20,Thioredoxin reductase,Prdx6
Accession No.	Uniprot ID:P99029
Uniprot	P99029
GeneID	54683;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prdx5

Note: This product is for in vitro research use only