

Pregnancy-specific beta-1-glycoprotein 11

Catalog No: #AP87000



Package Size: #AP87000-1 50ug #AP87000-2 100ug #AP87000-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pregnancy-specific beta-1-glycoprotein 11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PS-beta-G-11,Pregnancy-specific beta-1-glycoprotein 13,PS-beta-G-13,PSG13,PSG14
Accession No.	Uniprot ID:Q9UQ72
Uniprot	Q9UQ72
GeneID	5680;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PSG11

Note: This product is for in vitro research use only