

Protein patched homolog 2

Catalog No: #AP87001



Package Size: #AP87001-1 50ug #AP87001-2 100ug #AP87001-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein patched homolog 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PTC2,UNQ560/PRO1121/PRO57079
Accession No.	Uniprot ID:Q9Y6C5
Uniprot	Q9Y6C5
GeneID	8643;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PTCH2

Note: This product is for in vitro research use only