

Sulfhydryl oxidase 1

Catalog No: #AP87007



Package Size: #AP87007-1 50ug #AP87007-2 100ug #AP87007-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sulfhydryl oxidase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	hQSOX, Quiescin Q6, QSCN6, UNQ2520/PRO6013
Accession No.	Uniprot ID:O00391
Uniprot	O00391
GeneID	5768;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: QSOX1

Note: This product is for in vitro research use only