

Regenerating islet-derived protein 3-alpha

Catalog No: #AP87010



Package Size: #AP87010-1 50ug #AP87010-2 100ug #AP87010-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Regenerating islet-derived protein 3-alpha |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | REG-3-alpha, Islet of Langerhans regenerating protein 3, REG 3, Lithostathine 3, Pancreatitis-associated protein 2, RegIII, Regenerating islet-derived protein III-alpha, Reg III-alpha, Pap2, Reg3 |
| Accession No. | Uniprot ID:P35231 |
| Uniprot | P35231 |
| GeneID | 171162; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Reg3a

Note: This product is for in vitro research use only