

Serine/threonine-protein kinase PLK2

Catalog No: #AP87011



Package Size: #AP87011-1 50ug #AP87011-2 100ug #AP87011-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase PLK2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Polo-like kinase 2, PLK-2, Serine/threonine-protein kinase SNK, Serum-inducible kinase, Snk
Accession No.	Uniprot ID:Q9R012
Uniprot	Q9R012
GeneID	83722;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Plk2

Note: This product is for in vitro research use only