

Thyroid peroxidase

Catalog No: #AP87013



Package Size: #AP87013-1 50ug #AP87013-2 100ug #AP87013-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thyroid peroxidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	TPO
Accession No.	Uniprot ID:Q8HYB7
Uniprot	Q8HYB7
GeneID	403521;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TPO

Note: This product is for in vitro research use only