

Protein O-linked-mannose beta-1,2-N-acetylglucosaminyltransferase 1

Catalog No: #AP87016

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP87016-1 50ug #AP87016-2 100ug #AP87016-3 1mg

Description

Product Name	Protein O-linked-mannose beta-1,2-N-acetylglucosaminyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	POMGnT1
Accession No.	Uniprot ID:Q5XIN7
Uniprot	Q5XIN7
GeneID	362567;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pomgnt1

Note: This product is for in vitro research use only