

## Receptor-interacting serine/threonine-protein kinase 1

Catalog No: #AP87018



Package Size: #AP87018-1 50ug #AP87018-2 100ug #AP87018-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Receptor-interacting serine/threonine-protein kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cell death protein RIP, Receptor-interacting protein 1, RIP-1, Serine/threonine-protein kinase RIP, Rinp, Rip
Accession No.	Uniprot ID:Q60855
Uniprot	Q60855
GeneID	19766;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Ripk1

Note: This product is for in vitro research use only