

Pituitary homeobox 2

Catalog No: #AP87020



Package Size: #AP87020-1 50ug #AP87020-2 100ug #AP87020-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pituitary homeobox 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Homeobox protein PITX2,cPITX2,Paired-like homeodomain transcription factor 2,PTX2
Accession No.	Uniprot ID:O93385
Uniprot	O93385
GeneID	395862;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PITX2

Note: This product is for in vitro research use only