

Polycystin-1

Catalog No: #AP87021



Package Size: #AP87021-1 50ug #AP87021-2 100ug #AP87021-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Polycystin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Autosomal dominant polycystic kidney disease 1 protein homolog
Accession No.	Uniprot ID:O08852
Uniprot	O08852
GeneID	18763;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pkd1

Note: This product is for in vitro research use only