

## Phospholipid scramblase 1

Catalog No: #AP87027



Package Size: #AP87027-1 50ug #AP87027-2 100ug #AP87027-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phospholipid scramblase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PL scramblase 1,Ca(2+)-dependent phospholipid scramblase 1
Accession No.	Uniprot ID:P58195
Uniprot	P58195
GeneID	117540;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Plscr1

Note: This product is for in vitro research use only