

# Parathyroid hormone/parathyroid hormone-related peptide receptor

Catalog No: #AP87029

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP87029-1 50ug #AP87029-2 100ug #AP87029-3 1mg

## Description

Product Name	Parathyroid hormone/parathyroid hormone-related peptide receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	PTH/PTHrP type I receptor,PTH/PTHr receptor,Parathyroid hormone 1 receptor,PTH1 receptor,PTH1,PTHR1
Accession No.	Uniprot ID:Q9TU31
Uniprot	Q9TU31
GeneID	403779;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PTH1R

Note: This product is for in vitro research use only