

NMUR1 Antibody

Catalog No: #44375

Package Size: #44375-1 50ul #44375-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NMUR1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	NMUR1 Antibody detects endogenous levels of total NMUR1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human NMUR1
Target Name	NMUR1
Other Names	NMUR1, G protein-coupled receptor 66, GPR66, FM-3, FM3, G-protein coupled receptor 66, GPC-R, Neuromedin-U receptor 1, NMU-R1, NMU1R, Neuromedin U receptor 1
Accession No.	Swiss-Prot#: Q9HB89NCBI Gene ID:10316
Uniprot	Q9HB89
GeneID	10316;
Calculated MW	48kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

Western blotting: 1:1000

Images



Western Blot analysis of K562 cells using NMUR1 Polyclonal Antibody



Immunofluorescence analysis of LOVO cells, using NMUR1 Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only