

Proprotein convertase subtilisin/kexin type 9

Catalog No: #AP87035



Package Size: #AP87035-1 50ug #AP87035-2 100ug #AP87035-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proprotein convertase subtilisin/kexin type 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Neural apoptosis-regulated convertase 1,NARC-1,Proprotein convertase 9,PC9,Subtilisin/kexin-like protease PC9,Narc1
Accession No.	Uniprot ID:P59996
Uniprot	P59996
GeneID	298296;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pcsk9

Note: This product is for in vitro research use only