

Programmed cell death protein 2

Catalog No: #AP87036



Package Size: #AP87036-1 50ug #AP87036-2 100ug #AP87036-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Programmed cell death protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Zinc finger protein Rp-8,Rp8
Accession No.	Uniprot ID:P47816
Uniprot	P47816
GeneID	58934;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pdcd2

Note: This product is for in vitro research use only