

## Retinoblastoma-like protein 1

Catalog No: #AP87037



Package Size: #AP87037-1 50ug #AP87037-2 100ug #AP87037-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Retinoblastoma-like protein 1                          |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Hu   |
| Other Names        | 107 kDa retinoblastoma-associated protein,p107,pRb1    |
| Accession No.      | Uniprot ID:P28749                                      |
| Uniprot            | P28749   |
| GeneID             | 5933;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:RBL1

Note: This product is for in vitro research use only