

Renin receptor

Catalog No: #AP87039

Package Size: #AP87039-1 50ug #AP87039-2 100ug #AP87039-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Renin receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ATPase H(+)-transporting lysosomal accessory protein 2,ATPase H(+)-transporting lysosomal-interacting protein 2,ER-localized type I transmembrane adaptor,Embryonic liver differentiation factor 10,N14F,Renin/prorenin receptor,Vacuolar ATP synthase membrane sector-associated protein M8-9,ATP6M8-9,ATP6IP2,CAPER,ELDF10,HT028,MSTP009,PSEC0072
Accession No.	Uniprot ID:O75787
Uniprot	O75787
GeneID	10159;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATP6AP2

Note: This product is for in vitro research use only