

Platelet-derived growth factor C

Catalog No: #AP87043



Package Size: #AP87043-1 50ug #AP87043-2 100ug #AP87043-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet-derived growth factor C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PDGF-C, Fallotin, Spinal cord-derived growth factor, SCDGF, VEGF-E, Scdgc
Accession No.	Uniprot ID:Q8C119
Uniprot	Q8C119
GeneID	54635;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pdgfc

Note: This product is for in vitro research use only