

Platelet factor 4

Catalog No: #AP87044



Package Size: #AP87044-1 50ug #AP87044-2 100ug #AP87044-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Platelet factor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PF-4,C-X-C motif chemokine 4,CXCL4,SCYB4
Accession No.	Uniprot ID:P02777
Uniprot	P02777
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PF4

Note: This product is for in vitro research use only