

85/88 kDa calcium-independent phospholipase A2

Catalog No: #AP87045



Package Size: #AP87045-1 50ug #AP87045-2 100ug #AP87045-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	85/88 kDa calcium-independent phospholipase A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cal-PLA2,iPLA2,Group VI phospholipase A2,GVI PLA2
Accession No.	Uniprot ID:P97819
Uniprot	P97819
GeneID	53357;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pla2g6

Note: This product is for in vitro research use only