

Phosphotriesterase-related protein

Catalog No: #AP87047



Package Size: #AP87047-1 50ug #AP87047-2 100ug #AP87047-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphotriesterase-related protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Parathion hydrolase-related protein,hPHRP
Accession No.	Uniprot ID:Q96BW5
Uniprot	Q96BW5
GeneID	9317;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PTER

Note: This product is for in vitro research use only