

Ras-related protein Rab-14

Catalog No: #AP87049



Package Size: #AP87049-1 50ug #AP87049-2 100ug #AP87049-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ras-related protein Rab-14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:P61107
Uniprot	P61107
GeneID	94197;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rab14

Note: This product is for in vitro research use only