

Presenilin-2

Catalog No: #AP87052



Package Size: #AP87052-1 50ug #AP87052-2 100ug #AP87052-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Presenilin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PS-2,Ps2,Psnl2
Accession No.	Uniprot ID:O88777
Uniprot	O88777
GeneID	81751;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Psen2

Note: This product is for in vitro research use only