

Phosphatidylinositol phosphatase PTPRQ

Catalog No: #AP87053



Package Size: #AP87053-1 50ug #AP87053-2 100ug #AP87053-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphatidylinositol phosphatase PTPRQ
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Receptor-type tyrosine-protein phosphatase Q,PTP-RQ
Accession No.	Uniprot ID:Q9UMZ3
Uniprot	Q9UMZ3
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PTPRQ

Note: This product is for in vitro research use only